

DSM-TACE OF UNRESECTABLE MULTINODULAR BILOBAR HCC IN NAFLD CIRRHOSIS



Author: **Roberto Iezzi, MD**
Fondazione Policlinico Universitario A. Gemelli, IRCCS
Università Cattolica del Sacro Cuore, Rome, Italy



Patient

- 68 year old male
- Histologically confirmed **unresectable multinodular bilobar HCC G2-3** (>6 lesions, largest lesion 7 cm in hepatic segments IV-V) | Fig 1a-f
nonalcoholic fatty liver disease (NAFLD) cirrhosis
- BCLC B, Child-Pugh: A5 – MELD 7. ECOG 0. No ascites – no portal hypertension. Comorbidities: diabetes and arterial hypertension
- Lab parameters: **AFP** 5 ng/dl | **Total bilirubin** 0.5 mg/dl | **PLT** 195 x 10³/dl | **INR** 1.05 | **Albumin** 4.2g/dl | **Creatinine** 0.94 mg/dl
- Tumor board decision:
 - **DSM-TACE** with 50 mg Doxorubicin
 - **Bilobar treatment** (four treatments at 2-week interval; the first and third treatment were targeted to the lobe more involved by disease)




Figure 1: Pre-procedural CT scans show multinodular bilobar HCC lesions (up to 7 cm in size) (a-f)



DSM-TACE Procedure

- DSM-TACE procedures were performed in an angiographic suite, using patient monitoring and anesthesiological assistance under local anesthesia
- Selective lobar catheterization was performed with a coaxial technique placing a 2.7-Fr microcatheter in the right or left hepatic artery that was feeding the tumor lesions | Fig 2
- Solution of **450 mg in 7.5 ml of microspheres type EmboCept® S*** mixed with **50 mg of Doxorubicin** and non-ionic contrast medium was slowly infused in two steps:
 - Drug uptake: a ready-to-use solution composed of 50 mg of Doxorubicin diluted in 5 ml of saline solution plus 3.5 ml of EmboCept® S* plus 15 ml for right lobe or 10 ml for left lobe of non-ionic contrast medium was injected
 - Stop flow: the last 4 ml of EmboCept® S* plus 6 ml of non-ionic contrast medium was injected
- A technical success was obtained; in particular, all the expected dose was infused and a final “stop flow” was obtained




Figure 2: Angiogram shows lobar catheterization of right hepatic artery, using coaxial technique

DSM-TACE OF UNRESECTABLE MULTINODULAR BILOBAR HCC IN NAFLD CIRRHOSIS



Outcome

- Patient experienced no adverse events, with dismissal after 24 hours, without any pain or periprocedural complications
- 3-month CT follow-up showed an **almost complete response with necrosis of almost all multinodular hepatic lesions** | Fig 3a-f




Figure 3: 3-month post DSM-TACE CT scans show an almost complete response with necrosis of almost all hepatic lesions, without complications (a-f)



Outlook

- Based on 3-month follow-up result, patient will receive 2 more DSM-TACE sessions



CONCLUSION

- ▶ **DSM-TACE** with doxorubicin is a **safe, feasible, and effective option** in the treatment of patients with multinodular bilobar HCC, without drug-related toxicities
- ▶ Degradable Starch Microspheres cause a temporary occlusion, with a short ischemic period, allowing an **optimal drug uptake** with no post-embolic effects
- ▶ The transient vascular occlusion generated by **DSM allows to repeat treatment**, reducing the risk of liver toxicity that may occur when repeating cTACE or DEB-TACE
- ▶ DSM-TACE is used for **treatment of multinodular extensive intermediate HCC stage patients** where superselective TACE is not possible or could be characterized by a high risk of hepatic toxicities or adverse events. DSM-TACE is well-tolerated and effective for BCLC-B patients and BCLC-C patients ineligible for Sorafenib treatment

* Patient treated with EmboCep[†] S, which is equivalent to the successor and available product EmboCep[†] S DSM 50 µm [data on file].

DSM Degradable Starch Microspheres
TACE Transarterial chemoembolization

PharmaCept GmbH, Bessemerstr. 82, 12103 Berlin, Germany
Phone: +49-(0)30-7565985-0, Fax: +49-(0)30-7565985-11, info@pharmacept.com
EmboCep[®] is a registered trademark of PharmaCept GmbH.
© 2020 PharmaCept GmbH, Berlin, Germany.
202009004-09/2020



pharmacept.com